

Iron in low resolution? Removing molecular interferences with the *Neoma MS/MS MC-ICP-MS*.

GRANT CRAIG, MARKUS PFEIFER, JENNY ROBERTS,
CLAUDIA BOUMAN, NICHOLAS S. LLOYD AND
JOHANNES SCHWIETERS

Thermo Fisher Scientific (Bremen) GmbH

Presenting Author: grant.craig@thermofisher.com

Measuring the isotopic composition of iron, especially the isotope ratios $^{56}\text{Fe}/^{54}\text{Fe}$ and $^{57}\text{Fe}/^{54}\text{Fe}$, is a frequent and longstanding application for MC-ICP-MS instrumentation. The presence of common molecular interferences, including $^{40}\text{Ar}^{16}\text{O}^+$ and $^{40}\text{Ar}^{14}\text{N}^+$, has made high mass resolving power ($\geq 6,000 \Delta M/M$) a requirement to measure Fe isotope ratios [1]. However, high mass resolving power cuts ion transmission, increasing the sample size required to achieve high isotope ratio precision.

The re-introduction of MC-ICP-MS containing collision/reactions cells has opened the possibility of measuring Fe isotope ratios in low resolution. Wang *et al.* used H_2 and He gas in a collision/reaction cell to eliminate argon-based interferences on Fe and measure high precision $^{56}\text{Fe}/^{54}\text{Fe}$ in low resolution [2]. The Thermo Scientific™ Neoma MS/MS™ MC-ICP-MS includes the patented combination of a MC-ICP-MS with collision/reaction cell and pre-cell mass filter [3,4]. The groundbreaking and novel pre-cell mass filter, for fidelity of mass bias, has been described in Craig *et al.* [5]. Here we report experiments measuring $^{56}\text{Fe}/^{54}\text{Fe}$ and $^{57}\text{Fe}/^{54}\text{Fe}$ in low resolution using O_2 in the collision reaction cell. Using the pre-cell mass filter to prevent $^{40}\text{Ar}^+$ and $^{40}\text{Ar}^1\text{H}^+$ from entering the collision/reaction cell was necessary to allow O_2 to break up the Ar-based molecular interferences without adversely affecting Fe transmission: at the same time as enhancing Fe sensitivity. As a result, both $^{56}\text{Fe}/^{54}\text{Fe}$ and $^{57}\text{Fe}/^{54}\text{Fe}$ could be successfully measured free of molecular interferences in low resolution. The resulting method was applied to a selection of Fe samples, with long-term stability of $\delta^{56}\text{Fe}_{\text{IRMM-014}}$ and $\delta^{57}\text{Fe}_{\text{IRMM-014}}$ better than $\pm 0.02 \text{ ‰}$ (1SD) and 0.03 ‰ (1SD) respectively.

[1] Weyer, S.; Schwieters, J. B. (2003), *Int. J. Mass Spectrom.* 226, 355–368.

[2] Wang, J.; Tang, D.-M.; Su, B.-X.; Yuan, Q.-H.; Li, W.-J.; Gao, B.-Y.; Chen, K.-Y.; Bao, Z.-A.; Zhao, Y. (2022), *J. Anal. At. Spectrom.* 37, 1869–1875.

[3] Schwieters, J. B.; Elliott, T.; Coath, C. (2015), GB2546060B.

[4] Schwieters, J. B.; Jung, G. (2022) EP3769334B1.

[5] Craig, G.; Wehrs, H.; Bevan, D. G.; Pfeifer, M.; Lewis, J.; Coath, C. D.; Elliott, T.; Huang, C.; Lloyd, N. S.; Schwieters, J. B. (2021), *Anal. Chem.*, 93, 10519–10527.