

Chalcophile elements on collision course: S and Se isotope ratio measurements using the Sapphire Dual Path MC-ICP-MS with Collision/Reaction Cell

DANIEL PETERS AND YVAN GÉRARD

Nu Instruments Ltd

Presenting Author: daniel.peters@ametek.com

Chalcophile elements and their isotope ratios provide insights into a wide range of terrestrial processes from mantle fractionation to hydrothermal alteration and enable environmental and biogeochemical tracing. In this study, we will present isotope ratio data for S and Se measured on a Nu Sapphire Dual Path MC-ICP-MS with Collision/Reaction Cell to remove isobaric interferences.

Plasma-based S isotope measurements suffer from several isobaric interferences such as $^{16}\text{O}_2^+$ on $^{32}\text{S}^+$, $^{16}\text{O}^{17}\text{O}^+$ on $^{33}\text{S}^+$, $^{16}\text{O}^{18}\text{O}^+$ on $^{34}\text{S}^+$, and $^{18}\text{O}_2^+$ on $^{36}\text{S}^+$. For Se, numerous Ar-based (Ar dimers, $^{40}\text{Ar}^{37}\text{Cl}$) and Kr interferences affect all Se isotopes and have limited the exploration of this isotopic system thus far. While O, Ar, and Kr-based interferences can be well-suppressed using He-H₂ gas mixtures, both isotopic systems are prone to forming SH and SeH that will limit the accuracy of isotope ratio measurements. Thus, for both systems, a variety of reaction gases will be evaluated including He-H₂ and He-SF₆ gas mixtures as well as N₂.

Preliminary results using He-H₂ as reaction gas for Se reveal efficient Ar-Ar and Kr interference removal but significant SeH formation rates (~3%). After correction for SeH based on in-run measurements of $^{78}\text{SeH}^+$ relative to $^{78}\text{Se}^+$, Se isotope ratios in wet plasma using a self-aspirating nebuliser are well reproducible with a short-term repeatability of, e.g., $2\sigma_s$ SD < 35 ppm for $^{80}\text{Se}/^{78}\text{Se}$, $2\sigma_s$ SD < 45 ppm for $^{76}\text{Se}/^{78}\text{Se}$, or $2\sigma_s$ SD < 40 ppm for $^{77}\text{Se}/^{78}\text{Se}$ (with s = sample). Individual measurement errors are between 10-40 ppm ($2\sigma_s$ SE) for the different major Se isotope ratios. Isotope fractionation is thereby demonstrable mass-dependent with linear regression slopes of 1.994-1.997 ($R^2 = 0.999$) for $\ln(^{82}\text{Se}/^{78}\text{Se})$ vs. $\ln(^{80}\text{Se}/^{78}\text{Se})$ over two analytical sessions.