

Simultaneous pre-concentration of thorium and neodymium from seawater: Method development and application to the Kerguelen plateau

HABACUC PÉREZ-TRIBOULLIER^{1,3*}, Z. CHASE¹, T.L. NOBLE¹, A.T. TOWNSEND² AND A.R. BOWIE^{1,3}

¹Institute for Marine and Antarctic Studies, University of Tasmania, Hobart, Australia

*habacuc.pereztribouillier@utas.edu.au

²Central Science Laboratory, University of Tasmania, Hobart, Australia ³Antarctic Climate and Ecosystems Cooperative Research Centre, Hobart, TAS, Australia

Dissolved thorium (Th) and neodymium (Nd) are effective tracers of oceanic biogeochemical cycles and watermass mixing. Low concentrations of dissolved Th and Nd in seawater make their measurement a complex process, which requires pre-concentration and removal from the seawater matrix followed by further elemental purification, prior to detection using atomic spectroscopic methods. Here, we present a new method to simultaneously pre-concentrate Th and Nd using a chelating resin (NOBIAS) alongside modifications to existing separation techniques[1], with SF-ICP-MS detection. This method was applied to samples from the central Kerguelen plateau near to Heard and McDonald Islands, an area with strong natural iron fertilization and with a very dynamic hydrography.

Adsorption of trace elements by the NOBIAS resin is highly pH dependant[2]. We found the optimal pH to simultaneously pre-concentrate Th and Nd from 5-10 L of seawater was 4.5. At this pH the recovery of Th is > 75% and Nd > 90%. Our method showed good agreement with consensus values for GEOTRACES intercalibration materials (Th: BATS2000, SW2010-1; Nd: BATS2000, JNdi-1), with replicate measurements agreeing within 5%. Procedural blanks were less than 5 pg and 40 pg for Th and Nd respectively. Early results from the Kerguelen plateau indicate lower concentrations of ²³²Th than published data[3].

The accuracy, precision, blank levels and convenience of our method make it a good addition to existing techniques to measure Th and Nd in seawater, offering advantages of reduced processing time and sampling logistics.

[1] Anderson et al. (2012) *Limnol. Oceanogr. Methods* 10, 179-213. [2] Sohrin et al. (2008) *Anal. Chem.* 80(16), 6267-6273. [3] Venchiartutti et al. (2008) *Deep-Sea Res. I.* 55(10), 1343-1363.

**This abstract is too long to be accepted for publication.
Please revise it so that it fits into the column on one
page.**